

Date: 2011.05

Version 11_01

4X Lower Gel Buffer

Cat No : PT-4X-LOW-500ML Size : 500 mL

Description :

4X Lower Gel Buffer contain Tris and SDS at the appropriate pH and concentration for convenient and quick preparation of polyacrylamide gels for protein analysis by SDS-PAGE. Dilution of the Running Gel 4X buffer (1.5M Tris pH 8.8 and 0.4% SDS) to a 1X solution results in the standard, 0.375M Tris, 0.1% SDS concentrations routinely used in the preparation of the running gel for SDS-PAGE.

Component : *Store at RT*

4X Lower Gel Buffer (1.5M Tris pH 8.8 and 0.4% SDS)

Research Use Only

Please do not hesitate to contact us while you have any questions.

Manufactured for and distributed by Protech Technology Enterprise Co., Ltd.

TEL : +886-2-2655-7677 / FAX : +886-2-2655-7601 / Toll Free : 0800-231-530

E-MAIL: service@bio-protech.com.tw ; tech@bio-protech.com.tw

WEB: <http://www.bio-protech.com.tw>