

UltraPure DEPC Water

Catalog Number: PT-P560-500

Size : 500 mL

Storage : Room temperature

Description : UltraPure DEPC Water is a high quality reagent grade water produced by a professional water purification system. The water is treated with diethylpyrocarbonate (DEPC) and autoclaved to remove the residual DEPC. After production the UltraPure DEPC Water is dispensed into sterile RNase free containers and ready to use. It is tested for DNase, RNase and proteinase contamination through the standard method.

The UltraPure DEPC Water is indicated for reagent preparation in molecular biology, nucleic acid research, diagnostics, and for many other application requiring RNase free water.

Quality Control : DNase Activity ——— Non Detected

RNase Activity ——— Non Detected

Proteinase Activity ——— Non Detected

For Research Use Only.

Please do not hesitate to contact us if you have any questions.

Manufactured for and distributed by Protech Technology Enterprise Co.,Ltd
TEL : +886-2-2655-7677 / FAX : +886-2-2655-7601 / Toll Free : 0800-231-530
E-MAIL: service@bio-protech.com.tw ; tech@bio-protech.com.tw
WEB: <http://www.bio-protech.com.tw>