

UltraPure Water

Catalog Number: PT-P570

Size : 1 L

Storage : Room temperature

Description : UltraPure Water is a high quality reagent grade water produced by a professional water purification system. After production the water is sterile filtered and dispensed into sterile RNase free containers. The water is tested for DNase, RNase and proteinase contamination through the standard method.

The UltraPure Water is indicated for reagent preparation in molecular biology, nucleic acid research, diagnostics, tissue culture, *in vitro* fertilization, residue analysis, HPLC, spectrophotometry, and for many other application requiring exceptional quality water.

Quality Control : DNase Activity —— Non Detected

RNase Activity —— Non Detected

Proteinase Activity —— Non Detected

For Research Use Only.

Please do not hesitate to contact us if you have any questions.

Manufactured for and distributed by Protech Technology Enterprise Co.,Ltd

TEL : +886-2-2655-7677 / FAX : +886-2-2655-7601 / Toll Free : 0800-231-530

E-MAIL: service@bio-protech.com.tw ; tech@bio-protech.com.tw

WEB: <http://www.bio-protech.com.tw>

Protech Technology Enterprise Co., Ltd.